

# LUXTEC®



## ULTRALITE® PRO HEADLIGHT

Revolutionary. Engineered for Maximum Comfort and Stability.

 **INTEGRA™**  
Surgical

**THREE POINT WEIGHT ADJUSTMENT****DESIGNED FOR MAXIMUM COMFORT  
DURING EXTENDED WEAR****I OCCIPITAL BASKET ADJUSTMENT**

Allows optimum position of the occipital basket to balance the weight of the headlight and reduce neck strain.

**II CRANIAL PAD ADJUSTMENT**

Optimized location of cranial pad supports position of the head set and reduces temporal compression to a minimum.

**III HEAD BAND ADJUSTMENT**

Provides perfect fit and stability.



## IMPROVED LIGHT MODULE

### SUPERIOR OPTICS

- 10% increased light output\*
- Superior color fidelity

### ENHANCED DESIGN

- 10% weight reduction
- Horizontal movement eliminated
- Extended design for better vertical alignment with surgical loupes



## ULTRALITE® FIBER OPTIC CABLE

- Lightweight and flexible 3mm cable
- Bifurcated design allows even weight distribution and balance
- Translucent Sil-Kore™ tubing for durability and easy cleaning
- Provides visualization of fiber integrity



## MLX300 WATT XENON LIGHT SOURCE

- True white light gives high color fidelity
- New smart fan technology achieves quieter operation in the operating room
- Universal output for all fiberoptic cables
- **1000** hour warranty on Xenon bulb
- Proven quality and reliability



# THE INSTRUMENT PEOPLE

LUXTEC | OMNI-TRACT | JAIT | PADGETT | RB | REDMOND

## ORDERING INFORMATION

Part No.	Description
AX2000BIF	UltraLite® Pro headlight assembly with 9 ft (275cm) bifurcated cable, Light module and gown clips
AX2000BIFSL	UltraLite® Pro headlight assembly with 9 ft (275cm) bifurcated cable, Light module and gown clips with short linkage
001370LX	UltraLite® Pro replacement headband
001380LX	Light module, black
001380LXSL	Light module, black, short linkage
001378LX9	UltraLite® Pro bifurcated fiberoptic cable, 9 ft (275cm)
001378LX11	UltraLite® Pro bifurcated fiberoptic cable, 11 ft (335cm)
001378LX7P5	UltraLite® Pro bifurcated fiberoptic cable, 7.5 ft (230cm)
001387LX	Replacement Pad kit, UltraLite® Pro headlight
00MLX	MLX 300W Xenon Lightsource, with 20' power cord
00MLXEU	MLX 300W Xenon Lightsource, with 6.1m EU power cord
00MLXAU	MLX 300W Xenon Lightsource, with 6.1m Australian power cord
00MLXUK	MLX 300W Xenon Lightsource, with 6.1m UK power cord
001337	Floorstand for MLX 300W Xenon Lightsource



[www.Integra-Surgical.com](http://www.Integra-Surgical.com)

The Integra wave logo and Integra Surgical are trademarks of Integra LifeSciences Corporation.

SIL-KORE is a trademark of W.L. Gore & Associates, Inc.

Luxtec and UltraLite are registered trademarks of Integra LifeSciences Corporation.

© 2009 Integra LifeSciences Corporation